

vatech

EzRay AirTM Plus

WEIGHT INNOVATION



WORLD'S FIRST CNT PORTABLE X-RAY
FASTER WORKFLOW




WEIGHT INNOVATION WITH CNT (CARBON NANO TECHNOLOGY)

The world's first dental application of Carbon Nano Technology.

The EzRay Air Plus is a lightweight portable X-ray device designed for easy handling and stable positioning for optimal image quality on your intra-oral X-rays.



DOUBLE SCATTER SHIELD DESIGN FOR OPERATOR SAFETY



LOWER COOL DOWN TIMES, FASTER WORKFLOW

Carbon Nano Technology reduces cooldown time by up to 75%, letting users take more X-rays with less waiting and more focus on patient care.

NO WARMUP TIME

With no warm-up needed, Carbon Nano Technology enables instant exposure for a faster, smoother workflow.

SPECIFICATIONS [EzRay Air Plus]

Focal Spot	0.4 mm
Tube Voltage [kV]	65 kV
Tube Current [mA]	2.5 mA
Exposure Time Range	0.05 ~ 1.00 sec.
Total Filtration	Min. 1.5 mm Al
Source to Skin Distance	200 mm
Maximum Duty Cycle	1 : 60
Optional Accessories	Back Scatter Shield, Carry Case

* The specifications are subject to change without prior notice.