

PRECISION SIMPLIFIED

Apex Locator



Apex Locator

valt - Apex Locator


Simple & Intuitive design

Features :


- Up to 98% Accuracy
- Innovative resonance technique
- No calibration, No Zero setting (Quick to use)
- First in class Probe detection feature
- Stable and responsive graph
- Simplified display with sound feedback
(Sound frequency changes as you approach apical foramen with separate tone when you reach apex)
- Rechargeable Li-ion battery with long battery life (no need for frequent battery replacement)
- Auto power-off feature to enhance battery life
- Autoclavable Accessories
- Antimicrobial coating over the lip hook
- 1 year warranty on the main unit & six months warranty on accessories.

Standard product packaging includes

Innovopex- 1 Apex Locator Unit	- Qty 1
Power Adaptor & Cable	- Qty 1
Main cable	- Qty 1
Probe cable	- Qty 2
Lip Hook	- Qty 3
User Manual Warranty Card	



Measures root canals accurately and precisely



Technical Specifications :

Weight	150 GM.
Dimensions	100 X 90 X 90mm.
Power Consumption	< 0.3 W
Battery	3.7 V 1000 mAh, Li-Ion rechargeable battery
Charger Specs	5V, 2.4 amp, C Type