

vatech

WORLD'S No. 1 MOST RELIABLE

LARGEST SELLING

Intra Oral Sensor for SMART Digital Dentistry



ALSO AVAILABLE WITH
5 & 10 YEARS WARRANTY




Ez Sensor Classic

A Classic Choice of Classy Clinicians



Unique Size 1.5
Same Unit Suitable For
Adult & Paedo Patients



Cut Corners & Tapered Edges
Easy Positioning & Optimal
Patient Comfort, No Pain



Thinnest Sensor 4.8
CMOS APS Detector
easy to position in patient's mouth



Water & Dust Resistant
IP68 - Protects From Water
& Resists Dust



Fast & Efficient Workflow
Direct USB Interface, allows
fast & easy transfer of images



250+ Animated Videos
Unique Consultation Tool, Increases
Treatment Acceptance Ratio

ONLY COMPANY HAVING
DIRECT MANUFACTURING
& DIRECT SERVICE NETWORK



40,000+

SUCCESSFUL INSTALLATIONS IN INDIA

STILL COUNTING...

Superior Measurable Accurate Reliable Technology

SMART



EzSensor Classic

A Classic Choice of Classy Clinicians

Never Compromise on Your Reputation - Be Smart - Buy Branded

Redefining Digital Dentistry

**MOST RELIABLE BRAND OFFERING
VARIOUS WARRANTY OPTIONS**



Specification

Parameter	Size 1.5
Dimension [W x L x T]	39.5 x 29.2 x 4.8 mm
Active Area [W x L]	33.00 x 23.98 mm
Pixel Size	29.6 µm
Practical Resolution	17 lp/mm
Protection Code	IP 68, Dust & Water Proof
USB Cable length	2.7 meter

Complete Inhouse Manufacturing


HARDWARE & SOFTWARE

EzDent-i

Powerful Consultation Tool


- 2 in 1 communication with **EzDent-i**
Now it's possible to diagnose and consult at the same time on the same screen.
- **Unique 250+ Animated Consultation Videos by Vatech**
Quick Explanation of Procedure Increases Treatment Acceptance Ratio

Effective Diagnosis + Consultation



VIN ESC 0924

INTEGRATED
PMS



CleverDent

Patient Management Software


- **Fast Installation Setup** - Most Simple Design
- **Increased Productivity** - Faster Patient Consultation
- **Enhanced Reputation in Dentistry** - Get More Referrals
- **Easy Access of Data** - From Any Vatech Imaging Software

Widest Sales + SVC Network

450 + DIRECT PARTNERS ALL INDIA

vatech **DIRECT**
123 SVC
FASTEST SUPPORT MECHANISM

Register Your Concern on
www.vatechindia.in/direct



**TEAMVIEWER SUPPORT
FOR INSTANT SOLUTION**



Scan to Visit
Vatech Direct 123 SVC Page

VATECH INDIA PVT. LTD.

1016/2, 2nd Floor, S.M. Plaza, Bypass Road, Mahipalpur, New Delhi - 110037
W: www.vatechindia.in E: contactus@vatechindia.in M: 959 981 2206

#vatechindia

